

## **Diamond Polishing Paste- Instructions for Use**

### **Indications for use**

- Any products requiring a polished finish, including PMMA, acrylics, zirconia, porcelain, ceramics, and all metals

### **Composition**

- High diamond concentrate

### **How to Use**

- For a perfect surface texture, apply 1-2g of diamond polishing paste to each case using a polishing brush
  - A goat hair (19-21mm) or bison hair (18mm) polishing brush is recommended
- Recommended speed is 10,000-15,000 RPMs
- RG Diamond polish paste distributes evenly and smoothly, is splatter-proof, and won't leave behind a greasy residue

***Questions & Reordering: 855-472-8282***